

SET-UP

POWER SWITCHES (behind the table)	ON
KEITHLEY METER	PRESS DCI
HEATER (TTi)	POWER ON
HEATER (TTi) ON/OFF SWITCH	OFF
COMPUTER	TURN ON
PRESSURE GAUGE	POWER ON
SAMPLE	ALIGN
PRESSURE	APPLY
SAMPLE	(RE)ALIGN
CHILLER	WATER LEVEL AT LEDGE
CHILLER	POWER ON
CHILLER	PRESS 
CHILLER	WATCH FOR LEAKS
HEATER (TTi) ON/OFF SWITCH	ON

BEFORE TAKING IMAGE

NUC

PRESS

FOCUS

USE CLIP

TAKING IMAGE

RECORD

SAMPLE ID

TYPE OF LENS

BLIND

ALMOST CLOSE

WAIT

30 MINUTES

NUC

PRESS

RECORD

AMBIENT TEMPERATURE

HUMIDITY

PRESSURE

FILTER

AVERAGIN (GAIN 16)

TAKE IMAGE

RECORD

IMAGE NAME

REGION OF INTEREST

EXPORT

SHUT DOWN

HEATER (TTi) ON/OFF SWITCH	OFF
HEATER (TTi) POWER	OFF
CHILLER	OFF+UNPLUG
SAMPLE	UNLOAD
DATA	BACKUP
COMPUTER	SHUT DOWN
MAIN SIWTCHES (behind the table)	OFF